

# BIOCLEAN ANTI-FOG SPRAY For Autoclavable Goggles BAFS

## About

To address any fogging issues with our autoclavable goggles an application of BioClean Anti-fog Goggle Spray is recommended. The procedure for the successful application of the spray is as follows:

1. Autoclave the BioClean autoclavable goggles as usual
2. Just before donning, apply one or two sprays to the INSIDE of the lens
3. Using a sterile cleanroom wipe, dry the excess liquid from the inside surface of the lens with a spreading motion, covering the whole of the inside of the lens
4. Ensure the treated side of the lens is dry after wiping
5. Repeat the process after every autoclave cycle

## Specifications

**MATERIAL:**

Non-hazardous ingredients

**LENGTH:**

60ml

**COLOUR:**

Colourless

**STERILIZATION:**

Gamma Irradiated

## Features

- Reduces fogging
- Airless pump action spray
- Alcohol free
- Quick and easy application
- Suitable for all BioClean autoclavable goggles

ALCOHOL FREE

QUICK AND EASY APPLICATION





## Quality Standards

---

- Manufactured in a facility holding ISO 9001:2008 certifications



## Sterilization

---

Method: Gamma irradiation  
Minimum Dose: 25kGy

## Shelf Life & Storage

---

Three (3) years from date of manufacture. Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.

## TO ORDER

### Re-order code

BAFS-601

### Size

60ml

---

### Packing

One bottle per sealed inner PE bag;10 inner bags per sealed outer PE bag;10 outer bags per lined carton (100 bottles)

For additional information visit us at [www.bioclean.com](http://www.bioclean.com), or call us at +60 3 5569 3857/3859

### Nitritex (M) Sdn.Bhd

2 Jalan Jurunilai U1/20  
Seksyen U1, Hicom Glenmarie Industrial Park  
40150 Shah Alam, Selangor  
Malaysia

BAFS/1808/PDS11/APAC



Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of accessory for a specific application.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2018 Ansell Limited. All Rights Reserved.