

BIOCLEAN

ANTI-FOG SPRAY

For Autoclavable Goggles

BAFS

About

To address any fogging issues with our autoclavable goggles an application of BioClean Anti-fog Goggle Spray is recommended. The procedure for the successful application of the spray is as follows:

1. Autoclave the BioClean autoclavable goggles as usual
2. Just before donning, apply one or two sprays to the INSIDE of the lens
3. Using a sterile cleanroom wipe, dry the excess liquid from the inside surface of the lens with a spreading motion, covering the whole of the inside of the lens
4. Ensure the treated side of the lens is dry after wiping
5. Repeat the process after every autoclave cycle

Specifications

MATERIAL:

Non-hazardous ingredients

LENGTH:

60ml

COLOUR:

Colourless

STERILIZATION:

Gamma Irradiated

Features

- Reduces fogging
- Airless pump action spray
- Alcohol free
- Quick and easy application
- Suitable for all BioClean autoclavable goggles

ALCOHOL FREE

QUICK AND EASY APPLICATION



Billy Bioclean©



Quality Standards

- Manufactured in a facility holding ISO 9001:2008 certifications

Sterilization

Method: Gamma irradiation
Minimum Dose: 25kGy

Shelf Life & Storage

Three (3) years from date of manufacture. Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.

TO ORDER

Re-order code

BAFS-601

Size

60ml

Packing

One bottle per sealed inner PE bag; 10 inner bags per sealed outer PE bag; 10 outer bags per lined carton (100 bottles)

For additional information visit us at www.bioclean.com, or call us at +60 3 5569 3857/3859

Nitritex (M) Sdn.Bhd

2 Jalan Jurunilai U1/20
Seksyen U1, Hicom Glenmarie Industrial Park
40150 Shah Alam, Selangor
Malaysia

BAFS/1808/PDS11/APAC



Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of accessory for a specific application.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2018 Ansell Limited. All Rights Reserved.