

# NEW



# HAZ DETECT™

## HAZARDOUS DRUG MANIPULATION TECHNIQUE

TRAINING KIT

*"The OSHA HCS and USP chapter <800> require employee training for the tasks that will be performed as part of the (hazardous drug compounding) safety program."*

ASHP Guidelines on Handling Hazardous Drugs

Assists in training for drug handling competency and verification. Can be used to simulate spills and subsequent clean-up protocol processes.

Design a training method to reflect your most risk-associated process.

After performing the determined manipulation with the HAZ DETECT components, a black light shined on the working surface will reveal any spills and the need for technician retraining.

Cat. no. HVHAZ

### Kit Contents

- ▶ 5 x 50mg fluorescein powder septum vials
- ▶ 5 x 3ml reconstituted fluorescein septum tubes
- ▶ 5 tests per kit



*Needles, syringes, mini bag (see HVB10), and UV Lamp (see UVL56) not included.*