

CaviWipes™ HP

Surface Disinfectant

CaviWipes HP disinfecting wipes are the newest members of the trusted CaviWipes family providing a 1-minute universal contact time from a hydrogen peroxide formulation that qualifies for the EPA's rigorous Emerging Viral Pathogen Claim². The High Performance disinfection protocol offered by CaviWipes HP helps protect both patients and staff.



Efficient Disinfection

1-Step cleaning/disinfection³ and 1-minute universal contact time for all listed pathogens ensures efficient disinfection.



Broad Spectrum Efficacy

In addition to providing effective disinfection against all 67 listed pathogens including MDRO's (Multidrug Resistant Organisms), CaviWipes HP disinfecting wipes qualify for the EPA's rigorous Emerging Viral Pathogen claim², protecting you against future pathogens as well.



Ease of Use

The new green top opens and closes more easily and securely, and dispenses wipes with precision, helping to reduce waste.



Lowest Toxicity Level

US EPA Toxicity Category IV products minimize exposure to harmful cleaning chemicals. The use of gloves is not required.⁵



1 Minute Efficacy Against:⁴

Mycobacterium

Mycobacterium tuberculosis var: *bovis* (BCG)

Bacteria

Pseudomonas aeruginosa

Staphylococcus aureus

Salmonella enterica

Bordetella pertussis

Enterococcus faecium

Klebsiella aerogenes

Klebsiella pneumoniae

Escherichia coli

Escherichia coli O157:H7

Burkholderia cepacia

Stenotrophomonas maltophilia

Streptococcus pyogenes

Listeria monocytogenes

Shigella dysenteriae

Serratia marcescens

Bordetella bronchiseptica

Campylobacter jejuni

Proteus vulgaris

Multi-Drug-Resistant Bacteria (MDRO)

Penicillin-Resistant *Streptococcus pneumoniae* (PRSP)

Vancomycin-Resistant *Enterococcus faecalis* (VRE)

Methicillin-Resistant *Staphylococcus aureus* (MRSA)

Hospital-Acquired and Community Acquired Methicillin-Resistant *Staphylococcus aureus* (HA-MRSA and CA-MRSA)

Vancomycin Intermediate *Staphylococcus aureus* (VISA)

Vancomycin Resistant *Staphylococcus aureus* (VRSA)

Linezolid-Resistant *Staphylococcus aureus* (LRSA)

Carbapenem-Resistant *Klebsiella pneumoniae*

Multidrug-Resistant *Acinetobacter baumannii* (MRAB) Strains

ATCC 19606 and ATCC BAA-1605

NDM-1 *Enterobacter cloacae*

NDM-1 *Klebsiella pneumoniae*

MBL *Pseudomonas aeruginosa*

NMD-1 *Escherichia coli*

ESBL *Klebsiella pneumoniae*

ESBL *Escherichia coli*

Multidrug-Resistant *Klebsiella pneumoniae*

Carbapenem-Resistant *Escherichia coli*

Pathogenic Fungi

Trichophyton interdigitale

Candida auris

Candida albicans

Viruses (Non-Enveloped)

Human Norovirus

Human Rotavirus

Rhinovirus Type 37

Rhinovirus Type 1A

Human Adenovirus Type 5

Echovirus Type 11

Enterovirus Type 71

Viruses (Enveloped)

Herpes simplex virus Type 1 & Type 2

Avian Influenza A virus (H7N9)

Avian Influenza A virus (H5N1)

Influenza A virus (H3N2)

Influenza B virus (H1N1)

Parainfluenza virus Type 3

Hepatitis B virus (HBV)

Human Hepatitis C virus (HCV)

Human Immunodeficiency virus Type 1 (HIV-1)

Severe Acute Respiratory Syndrome-related

Coronavirus 2 (SARS-CoV-2)

Vaccinia virus

SARS-Associated Coronavirus (SARs-CoV)

Cytomegalovirus (CMV)

Influenza A virus (H1N1) Strains A/4108/2009 and A/04/09

Human Coronavirus

Respiratory Syncytial Virus (RSV)

Measles virus

Active Ingredient

Hydrogen Peroxide 0.51%

Other Ingredients† 99.49%

Total 100.00%

[† Does not include the weight of the towelette]



Ordering Information

Part No. Description

16-1100 CaviWipes HP (6" x 6.75") – 160 towelettes per canister

16-1150 CaviWipes HP XL (9" x 12") – 65 towelettes per canister

- List N: Disinfectants for Use Against SARS-CoV-2. Feb 28, 2022, www.epa.gov/pesticide-registration/list-n-disinfectants-use-against-sars-cov-2.
- CaviWipes HP qualifies for the EPA Emerging Viral Pathogen program against all virus tiers (enveloped, small and large non-enveloped).
- Follow Product Label for special precleaning and PPE instructions.
- Always follow Product Label for Cleaning/ Disinfection Directions.

To learn more about CaviWipes HP scan the QR code or visit metrex.com

